

Anti-GFP/eGFP-Tag Rabbit Monoclonal Antibody

Catalog Number: M30939-3


Overview

Product Name	Anti-GFP/eGFP-Tag Rabbit Monoclonal Antibody
Reactive Species	GFP Tag
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of GFP/eGFP-Tag protein. Validated for research with WB,IF,IP,ELISA applications. Anti GFP/eGFP-Tag is reactive for GFP-Tag samples
Application	ELISA, IP, IF, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P42212

Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	WB 1:5000-1:20000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200

Anti-GFP/eGFP-Tag Rabbit Monoclonal Antibody (M30939-3) Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-GFP/eGFP antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: 293F transfection by EGFP Predicted band size: 27kDa Observed band size: 27kDa

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GFP/eGFP-Tag Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.