

Anti-ECRG4 Mouse Monoclonal Antibody [Clone ID: OTI3B10]

Catalog Number: M32743


Overview

Product Name	Anti-ECRG4 Mouse Monoclonal Antibody [Clone ID: OTI3B10]
Reactive Species	Human
Description	Boster Bio C2orf40 mouse monoclonal antibody, clone OTI3B10 (formerly 3B10). Catalog# M32743. Tested in FC, IF, WB. This antibody reacts with Human.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI3B10
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9H1Z8


Technical Details

Immunogen	Full length human recombinant protein of human C2orf40 (NP_115787) produced in HEK293T cell.
Isotype	IgG2b
Concentration	0.85 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IF 1:100 Flow Cytometry 1:100


Anti-ECRG4 Mouse Monoclonal Antibody [Clone ID: OTI3B10] (M32743) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY C2orf40 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-C2orf40.



Flow cytometric Analysis of HeLa cells



Flow cytometric Analysis of Jurkat cells



Immunofluorescent staining of HeLa cells using anti-C2orf40 mouse monoclonal antibody (M32743).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ECRG4 Mouse Monoclonal Antibody [Clone ID: OTI3B10]

For Research Use Only. Not for use in diagnostic procedures.