

Anti-Spike RBD Reference Antibody (Regdanvimab)

Catalog Number: M34012-3

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Spike RBD Reference Antibody (Regdanvimab)
Reactive Species	Influenza A virus
Description	Boster Bio Anti-Spike RBD Reference Antibody (Regdanvimab) (Catalog # M34012-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Influenza A virus. Endotoxin: < 0.1014EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.


Technical Details

Isotype	IgG1
Form	Liquid


Anti-Spike RBD Reference Antibody (Regdanvimab) (M34012-3) Images



Anti-Spike RBD Reference Antibody (Regdanvimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-Spike RBD Reference Antibody (Regdanvimab) is more than 100%



Immobilized 2019nCoV RBD His at 2 μg/mL can bind Anti-Spike RBD Reference Antibody (Regdanvimab)

Anti-Spike RBD Reference Antibody (Regdanvimab)

For Research Use Only. Not for use in diagnostic procedures.