

Anti-TCR Reference Antibody (TOL101)

Catalog Number: M34015


Overview

Product Name	Anti-TCR Reference Antibody (TOL101)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TCR Reference Antibody (TOL101) (Catalog # M34015). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Human, Mouse. Endotoxin: < 0.1021EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.


Technical Details

Isotype	IgG2SA
Form	Liquid

Anti-TCR Reference Antibody (TOL101) (M34015) Images



Anti-TCR Reference Antibody (TOL101) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TCR Reference Antibody (TOL101) is more than 98.12%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TCR Reference Antibody (TOL101)

For Research Use Only. Not for use in diagnostic procedures.