

Anti-PCLA Reference Antibody (U.Penn. patent anti-PCLA)

Catalog Number: M34037

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-PCLA Reference Antibody (U.Penn. patent anti-PCLA)
Reactive Species	Human
Description	Boster Bio Anti-PCLA Reference Antibody (U.Penn. patent anti-PCLA) (Catalog # M34037). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.1068EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q9HGS6

Technical Details

Isotype	IgG1
Form	Liquid

Anti-PCLA Reference Antibody (U.Penn. patent anti-PCLA) (M34037) Images



Anti-PCLA Reference Antibody (U.Penn. patent anti-PCLA) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%

[Anti-PCLA Reference Antibody \(U.Penn. patent anti-PCLA\)](#)

For Research Use Only. Not for use in diagnostic procedures.