

## Anti-Tenascin TNN Antibody (Monoclonal, BC-24)

Catalog Number: MA1094

### About TNN

The tenascins are a family of extracellular matrix proteins (ECMs). Tenascin is an extracellular matrix protein. It expressed in an unusually restricted pattern during embryogenesis and has been implicated in a variety of morphogenetic phenomena. Human hexabrachion gene (tenascin) is mapped to chromosome 9, bands q32-q34. The coding region of the hexabrachion gene spans approximately 80 kilobases of DNA and consists of 27 exons separated by 26 introns.

### Overview

Product Name	Anti-Tenascin TNN Antibody (Monoclonal, BC-24)
Reactive Species	Human
Description	Boster Bio Anti-Tenascin TNN Antibody (Monoclonal, BC-24) catalog # MA1094. Tested in IHC, ICC, WB applications. This antibody reacts with Human.
Application	IHC, ICC, WB
Clonality	Monoclonal BC-24
Formulation	Mouse IgG in stabilizing components, 1.2% sodium acetate and 0.01mg NaN3.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UQP3

### Technical Details


Immunogen	Human tenascin.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Mouse IgG (EK1001) for Western blot, and HRP Conjugated anti-Mouse IgG Super Vision Assay Kit (SV0001-1) for IHC(P), IHC(F) and ICC.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Mouse IgG1
Form	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 ug/ml.
Purification	Ascites

Suggested Dilutions

Immunohistochemistry (Paraffin-embedded Section), 2-4ug/ml, Human  
Immunocytochemistry , 1ug/ml, Human, -  
Immunohistochemistry (Frozen Section), 2-4ug/ml, Human, -  
Western blot, 1-2ug/ml, Human

## Anti-Tenascin TNN Antibody (Monoclonal, BC-24) (MA1094) Images

---



Anti-Tenascin antibody (monoclonal), MA1094, IHC(P)IHC(P):  
Human Mammary Cancer Tissue

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tenascin TNN Antibody (Monoclonal, BC-24)

For Research Use Only. Not for use in diagnostic procedures.