

Anti-Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Monoclonal Antibody

Catalog Number: MP01723


Overview

Product Name	Anti-Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Monoclonal Antibody catalog # MP01723. Tested in WB, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal ADF-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O96013/Q9P286/Q9NQ5

Technical Details

Immunogen	A synthesized peptide derived from human Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, growth, proliferation or cell survival.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 FC 1:200

Anti-Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Monoclonal Antibody (MP01723) Images



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) expression in U87-MG cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.