

Anti-Firefly Luciferase Rabbit Monoclonal Antibody

Catalog Number: MT0021


Overview

Product Name	Anti-Firefly Luciferase Rabbit Monoclonal Antibody
Reactive Species	Firefly
Description	Boster Bio Anti-Firefly Luciferase Rabbit Monoclonal Antibody catalog # MT0021. Tested in WB, ICC/IF applications. This antibody reacts with Firefly.
Application	IF, ICC, WB
Clonality	Monoclonal DAB-6
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08659


Technical Details

Immunogen	A synthesized peptide derived from Firefly Luciferase
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:1000-1:5000 ICC/IF 1:50-1:200

Anti-Firefly Luciferase Rabbit Monoclonal Antibody (MT0021) Images



Western blot analysis of Firefly Luciferase Antibody expression in (1) Firefly Luciferase transfected 293T cell lysate; (2) HeLa cell lysate; (3) NIH/3T3 cell lysate; (4) C6 cell lysate.



Immunofluorescent analysis of 293 cells transfected with Firefly Luciferase, using Firefly Luciferase Antibody .

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Firefly Luciferase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.