

## Anti-p53 (acetyl K382) TP53 Rabbit Monoclonal Antibody

Catalog Number: P00001-2


### Overview

Product Name	Anti-p53 (acetyl K382) TP53 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-p53 (acetyl K382) TP53 Rabbit Monoclonal Antibody catalog # P00001-2. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Monoclonal EAE-20
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04637

### Technical Details

Immunogen	A synthesized peptide derived from human p53 (acetyl K382)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200

## Anti-p53 (acetyl K382) TP53 Rabbit Monoclonal Antibody (P00001-2) Images



Western blot analysis of p53 (acetyl K382) expression in Jurkat cell lysate treated with etoposide and TSA.

## 9 Publications Citing This Product

1. PubMed ID: 25395712, Li W, Wu D, Wei B, Wang S, Sun H, Li X, Zhang F, Zhang C, Xin Y. Afr J Tradit Complement Altern Med. 2014 Aug 23;11(5):99-104. Ecollection 2014. Anti-Tumor Effect Of Cactus Polysaccharides On Lung Squamous Carcinoma Cells (Sk-Mes-1).
2. PubMed ID: 26137081, Lactotransferrin expression is downregulated and affects the mitogen-activated protein kinase pathway in gastric cancer
3. PubMed ID: 27065079, TALENs-directed?knockout?of?the?full-length?transcription?factor?Nrf1?that?represses?malignantbehaviour?of?human?hepatocellular?carcinoma?(HepG2)?cells

Visit [bosterbio.com/anti-p53-acetyl-k382-rabbit-monoclonal-antibody-p00001-2-boster.html](http://bosterbio.com/anti-p53-acetyl-k382-rabbit-monoclonal-antibody-p00001-2-boster.html) to see all 9 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p53 (acetyl K382) TP53 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.