

## Anti-Phospho-EGFR (Y1068) Antibody

Catalog Number: P00023-7

### About EGFR

The epidermal growth factor receptor (EGFR; ErbB-1; HER1 in humans) is the cell-surface receptor for members of the epidermal growth factor family (EGF-family) of extracellular protein ligands. It is a member of the ErbB family of receptors, a subfamily of four closely related receptor tyrosine kinases: EGFR (ErbB-1), HER2/c-neu (ErbB-2), Her 3 (ErbB-3) and Her 4 (ErbB-4). EGFR exists on the cell surface and is activated by binding of its specific ligands, including epidermal growth factor and transforming growth factor alpha (TGFalpha). EGFR and its ligands are cell signaling molecules involved in diverse cellular functions, including cell proliferation, differentiation, motility, and survival, and in tissue development. Mutations that lead to EGFR overexpression (known as upregulation) or overactivity have been associated with a number of cancers, including lung cancer and glioblastoma multiforme. In this latter case a more or less specific mutation of EGFR, called EGFRvIII is often observed.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Phospho-EGFR (Y1068) Antibody                                                                                                                                                                                                                                                                                                                                                                                         |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                          |
| Description          | Boster Bio Anti-Phospho-EGFR (Y1068) Antibody catalog # P00023-7. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                    |
| Application          | IF, IHC, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                           |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                             |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Uniprot ID           | P00533                                                                                                                                                                                                                                                                                                                                                                                                                     |

### Technical Details

|                     |                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------|
| Immunogen           | A synthesized peptide derived from human EGFR around the phosphorylation site of Tyr1097.                      |
| Form                | Liquid                                                                                                         |
| Concentration       | 500 ug/ml                                                                                                      |
| Purification        | Protein A affinity purified.                                                                                   |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-200<br>Immunocytochemistry/Immunofluorescence, 1:50-200 |

[Anti-Phospho-EGFR \(Y1068\) Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.