

Anti-Phospho-STAT1 (Y701) Rabbit Monoclonal Antibody

Catalog Number: P00036


Overview

Product Name	Anti-Phospho-STAT1 (Y701) Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Phospho-STAT1 (Y701) Rabbit Monoclonal Antibody catalog # P00036. Tested in WB application. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal HHB-19
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42224


Technical Details

Immunogen	A synthesized peptide derived from human STAT1
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

Anti-Phospho-STAT1 (Y701) Rabbit Monoclonal Antibody (P00036) Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Phospho-STAT1 (Y701) expression in (1) A431 cell lysate; (2) A431 lysate cell treated with EGF.

2 Publications Citing This Product

1. PubMed ID: 28903386, Synergistic inhibition of colon cancer growth by the combination of methylglyoxal and silencing of glyoxalase I mediated by the STAT1 pathway
2. PubMed ID: 17696225, Hepatitis C virus non-structural 5A abrogates signal transducer and activator of transcription-1 nuclear translocation induced by IFN-alpha through dephosphorylation.

Visit bosterbio.com/anti-phospho-stat1-y701-rabbit-monoclonal-antibody-p00036-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-STAT1 (Y701) Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.