

Anti-Phospho-ER/ESR1 (Y537) Antibody

Catalog Number: P00057-2

About ESR1

Estrogen receptor alpha (ER-alpha), also known as NR3A1, is one of two main types of estrogen receptor, a nuclear receptor that is activated by the sex hormone estrogen. Estrogen receptors are involved in pathological processes including breast cancer, endometrial cancer, and osteoporosis. In humans, ER-alpha is encoded by the gene ESR1 (Estrogen Receptor 1). It is mapped to 6q25.1. This gene is a ligand-activated transcription factor composed of several domains important for hormone binding, DNA binding, and activation of transcription. The protein localizes to the nucleus where it may form a homodimer or a heterodimer with estrogen receptor 2. Estrogen and its receptors are essential for sexual development and reproductive function, it also play a role in other tissues such as bone.

Overview

Product Name	Anti-Phospho-ER/ESR1 (Y537) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-ER/ESR1 (Y537) Antibody catalog # P00057-2. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P03372

Technical Details

Immunogen	A synthesized peptide derived from human Estrogen Receptor- alpha around the phosphorylation site of Tyr537.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
Immunohistochemistry, 1:50-200
Immunocytochemistry/Immunofluorescence, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-ER/ESR1 (Y537) Antibody

For Research Use Only. Not for use in diagnostic procedures.