

Anti-Phospho-NRF2/NFE2L2 (S344) Antibody

Catalog Number: P00078-2

About NFE2L2

NFE2L2(nuclear factor (erythroid-derived 2)-like 2) also known as NRF2 or NFE2-RELATED TRANSCRIPTION FACTOR 2, is a transcription factor that in humans is encoded by the NFE2L2 gene. NFE2, NFE2L1, and NFE2L2 comprise a family of human genes encoding basic leucine zipper (bZIP) transcription factors. NFE2L2 induces the expression of various genes including those that encode for several antioxidant enzymes, and it may play a physiological role in the regulation of oxidative stress. The NFE2L2 gene is located on 2q31.2. The identification of somatic mutations that disrupt the NRF2-KEAP1 interaction to stabilize NRF2 and increase the constitutive transcription of NRF2 target genes indicated that enhanced reactive oxygen species (ROS) detoxification and additional NRF2 functions may in fact be tumorigenic. Oncogene-directed increased expression of Nrf2 is a mechanism for the activation of the Nrf2 antioxidant program evident in primary cells and tissues of mice expressing KRas(G12D) and BRaf(V619E), and in human pancreatic cancer.

Overview

Product Name	Anti-Phospho-NRF2/NFE2L2 (S344) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-NRF2/NFE2L2 (S344) Antibody catalog # P00078-2. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q16236

Technical Details

Immunogen	A synthesized peptide derived from human nrf2 around the phosphorylation site of Ser344.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
Immunohistochemistry, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-NRF2/NFE2L2 (S344) Antibody

For Research Use Only. Not for use in diagnostic procedures.