

## Anti-Phospho-PDGF Receptor beta (Y740) PDGFRB Rabbit Monoclonal Antibody

Catalog Number: P00096

**Submit a product review to Biocompare.com**

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Phospho-PDGF Receptor beta (Y740) PDGFRB Rabbit Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                           |
| Reactive Species     | Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Description          | Boster Bio Anti-Phospho-PDGF Receptor beta (Y740) PDGFRB Rabbit Monoclonal Antibody catalog # P00096. Tested in WB, ICC/IF applications. This antibody reacts with Mouse, Rat.                                                                                                                                                                                                                                                                     |
| Application          | IF, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Clonality            | Monoclonal FDF-16                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | P09619                                                                                                                                                                                                                                                                                                                                                                                                                                             |

### Technical Details

|               |                                                                            |
|---------------|----------------------------------------------------------------------------|
| Immunogen     | A synthesized peptide derived from human Phospho-PDGF Receptor beta (Y740) |
| Isotype       | Rabbit IgG                                                                 |
| Form          | Liquid                                                                     |
| Concentration | 0.5mg/ml                                                                   |

|                     |                                  |
|---------------------|----------------------------------|
| Purification        | Affinity-chromatography          |
| Suggested Dilutions | WB 1:500-2000<br>ICC/IF 1:50-200 |

## Anti-Phospho-PDGF Receptor beta (Y740) PDGFRB Rabbit Monoclonal Antibody (P00096) Images



Western blot analysis of Phospho-PDGF Receptor beta (Y740) expression in NIH 3T3 treated with PDGF.

### 1 Publications Citing This Product

1. PubMed ID: 18466260, RNA interference targeting the platelet-derived growth factor receptor  $\beta$  subunit ameliorates experimental hepatic fibrosis in rats

Visit [bosterbio.com/anti-phospho-pdgf-receptor-beta-y740-rabbit-monoclonal-antibody-p00096-boster.html](http://bosterbio.com/anti-phospho-pdgf-receptor-beta-y740-rabbit-monoclonal-antibody-p00096-boster.html) to see all 1 publications.

[Anti-Phospho-PDGF Receptor beta \(Y740\) PDGFRB Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.