

Anti-Phospho-FAK (Y397) PTK2 Rabbit Monoclonal Antibody

Catalog Number: P00151


Overview

Product Name	Anti-Phospho-FAK (Y397) PTK2 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-FAK (Y397) PTK2 Rabbit Monoclonal Antibody catalog # P00151. Tested in WB, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, ICC, WB
Clonality	Monoclonal EFO-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q05397

Technical Details

Immunogen	A synthesized peptide derived from human Phospho-FAK (Y397)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200

Anti-Phospho-FAK (Y397) PTK2 Rabbit Monoclonal Antibody (P00151) Images



Western blot analysis of Phospho-FAK (Y397) expression in EGF treated 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

4 Publications Citing This Product

1. PubMed ID: 10.3389/fphys.2018.01343, Knee Extensors Muscle Plasticity Over a 5-Years Rehabilitation Process After Open Knee Surgery
2. PubMed ID: 10.1016/j.yexmp.2019.03.006, Focal adhesion kinase coordinates costamere-related JNK signaling with muscle fiber transformation after Achilles tenotomy and tendon reconstruction
3. PubMed ID: 30337877, Knee extensors muscle plasticity over a 5-years rehabilitation process after open knee surgery

Visit bosterbio.com/anti-phospho-fak-y397-rabbit-monoclonal-antibody-p00151-boster.html to see all 4 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-FAK (Y397) PTK2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.