

Anti-Phospho-CHEK2 T68 CHEK2 Antibody

Catalog Number: P00277

About CHEK2

CHK2 (also known as CHEK2, Protein kinase CHK2 isoform a, and checkpoint-like protein) is a serine/ threonine-protein kinase involved in the control of cell cycle checkpoints and may also participate in transduction of the DNA damage and replicational stress signals. CHK2 is the mammalian ortholog of the budding yeast Rad53 and fission yeast Cds1 checkpoint kinases. The amino-terminal domain of CHK2 contains a series of seven serine and threonine residues (Ser19, Thr26, Ser28, Ser33, Ser35, Ser50 and Thr68) followed by glutamine (SQ or TQ motif). These are known to be preferred sites for phosphorylation by ATM/ATR kinases. Indeed, after DNA damage by ionizing radiation (IR), UV irradiation or hydroxyurea treatment, Thr68 and other sites in this region become phosphorylated by ATM/ATR. The SQ/TQ cluster domain, therefore, seems to have a regulatory function. Phosphorylation at Thr68 is a prerequisite for the subsequent activation step, which is attributable to autophosphorylation of Chk2 on residues Thr383 and Thr387 in the activation loop of the kinase domain. CHK2 inhibits CDC25C phosphatase by phosphorylating it on Ser-216, preventing the entry into mitosis. This kinase may have a role in meiosis as well. Kinase activity is up regulated by autophosphorylation and the protein is rapidly phosphorylated in response to DNA damage and to replication block.

Overview


Product Name	Anti-Phospho-CHEK2 T68 CHEK2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Phospho-CHEK2 T68 CHEK2 Antibody (Catalog # P00277). Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, IP, WB
Clonality	Polyclonal
Formulation	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.01% (w/v) Sodium Azide
Storage Instructions	Store vial at -20°C prior to opening. Aliquot contents and freeze at -20°C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of opening. (Ship on dry ice.)
Host	Rabbit
Uniprot ID	O96017

Technical Details

Immunogen	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to an internal region near aa 50-75 of Human CHK2.
Predicted Reactive Species	Chimpanzee

Isotype	IgG
Form	Liquid (sterile filtered)
Concentration	1.0mg/ml by UV absorbance at 280 nm
Purification	This affinity purified antibody is directed against the phosphorylated form of human CHK2 at the pT68 residue. The product was affinity purified from monospecific antiserum by immunoaffinity purification. Antiserum was first purified against the phosphorylated form of the immunizing peptide. The resultant affinity purified antibody was then cross adsorbed against the non-phosphorylated form of the immunizing peptide. Reactivity occurs against human CHK2 pT68 protein and the antibody is specific for the phosphorylated form of the protein. Reactivity with non-phosphorylated human CHK2 is minimal by ELISA. The antibody does not cross-react with Chk2 phosphorylated at other sites. A BLAST analysis was used to suggest reactivity with this protein from human and chimpanzee based on 100% homology for the immunogen sequence. Cross-reactivity with CHK2 protein from mouse and rat may occur as sequence homology varies by one amino acid residues in this sequence (90% homology). Cross-reactivity with CHK2 homologues from other sources has not been determined.
Suggested Dilutions	ELISA: 1:10,000 - 1:70,000 IP: 1:100 WB: 1:200 - 1:2,000 This affinity purified antibody has been tested for use in ELISA and by western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 61 kDa in size corresponding to CHK2 by western blotting in the appropriate cell lysate or extract. Less than 1% reactivity is observed against the non-phosphorylated form of the immunizing peptide. This antibody is phospho specific for pT68 of CHK2.

Anti-Phospho-CHK2 T68 CHEK2 Antibody (P00277) Images



Western blot using Boster's Affinity Purified anti-Chk2 pT68 antibody shows detection of a predominant band at ~60 kDa corresponding to phosphorylated Chk2 (arrowhead) in MCF-7 whole cell lysates after treatment with doxorubicin. Chk2 phosphorylation was induced using increasing concentrations of the DNA damaging agent doxorubicin as indicated for 24 h prior to lysate production. Personal communication, Xiao HeYang, University of Oklahoma Health Sciences Center.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-CHK2 T68 CHEK2 Antibody

For Research Use Only. Not for use in diagnostic procedures.