

## Anti-Phospho-PI3 Kinase p85 alpha (Y607) Antibody

Catalog Number: P00318

### Overview

Product Name	Anti-Phospho-PI3 Kinase p85 alpha (Y607) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-PI3 Kinase p85 alpha (Y607) Antibody catalog # P00318. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P27986

### Technical Details

Immunogen	A synthesized peptide derived from human PI3 Kinase p85 alpha around the phosphorylation site of Y607
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:1000 IHC 1:50-1:100

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-PI3 Kinase p85 alpha (Y607) Antibody

For Research Use Only. Not for use in diagnostic procedures.