

Anti-Phospho-NFATC1 (S172) Antibody

Catalog Number: P00340

About NFATC1

NFATC1(Nuclear factor of activated T-cells, cytoplasmic 1), also known as NFATC or NFAT2, is a protein that in humans is encoded by the NFATC1 gene. And it is also a component of the nuclear factor of activated T cells DNA-binding transcription complex. The NFATC1 gene is mapped on 18q23. Proteins belonging to this family of transcription factors play a central role in inducible gene transcription during immune response. NFATC1 was not detected in brain, liver, or kidney. NFATC1 is an inducible nuclear component which functions as a major molecular target for the immunosuppressive drugs such as cyclosporin A. The solution structure of the binary complex formed between the core DNA-binding domain of human NFATC1 and the 12-bp oligonucleotide duplex containing the ARRE2 DNA site from the IL2 promoter. Nfatc1 was expressed exclusively in mouse hair follicle stem cells, and using gain- and loss-of-function approaches, and Nfatc1 inhibited stem cell activation.

Overview

Product Name	Anti-Phospho-NFATC1 (S172) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-NFATC1 (S172) Antibody catalog # P00340. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O95644

Technical Details

Immunogen	A synthesized peptide derived from human NFAT2 around the phosphorylation site of Ser172.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
Immunohistochemistry, 1:50-200
Immunocytochemistry/Immunofluorescence, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-NFATC1 (S172) Antibody

For Research Use Only. Not for use in diagnostic procedures.