

## Anti-Phospho-LRP5 (Y1473) Antibody

Catalog Number: P00504

### About LRP5

Low-density lipoprotein receptor-related protein 5 is a protein that in humans is encoded by the LRP5 gene. It is mapped to 11q13.2. LRP5 is a transmembrane low-density lipoprotein receptor that binds and internalizes ligands in the process of receptor-mediated endocytosis. This protein also acts as a co-receptor with Frizzled protein family members for transducing signals by Wnt proteins and was originally cloned on the basis of its association with diabetes mellitus type 1 in humans. This protein plays a key role in skeletal homeostasis. In addition to that, the binding of axin to LRP5 is an important part of the Wnt signal transduction pathway, and it also acts as a target for the inhibitory effects of Dickkopf, another developmental protein, on Wnt signaling.

### Overview

Product Name	Anti-Phospho-LRP5 (Y1473) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-LRP5 (Y1473) Antibody catalog # P00504. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O75197

### Technical Details

Immunogen	A synthesized peptide derived from human LRP5 around the phosphorylation site of Tyr1473.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Immunohistochemistry, 1:50-200

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-LRP5 (Y1473) Antibody

For Research Use Only. Not for use in diagnostic procedures.