

Anti-Phospho-TGFR1 (T204) Antibody

Catalog Number: P00832-1

About TGFR1

Transforming growth factor, beta receptor I is a TGF beta receptor. TGFR1 is its human gene. The protein encoded by this gene forms a heteromeric complex with type II TGF-beta receptors when bound to TGF-beta, transducing the TGF-beta signal from the cell surface to the cytoplasm. Mutations in this gene have been associated with Loeys-Dietz aortic aneurysm syndrome (LDAS). TGFR1 regulates cell cycle progression by a unique signaling mechanism that involves its binding to TGFR2 and activation of TGFR1. Both are transmembrane serine/threonine receptor kinases. The TGFR1 receptor may be inactivated in many of the cases of human tumor cells refractory to TGF-beta-mediated cell cycle arrest.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Phospho-TGFR1 (T204) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-TGFR1 (T204) Antibody catalog # P00832-1. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P36897

Technical Details

Immunogen	A synthesized peptide derived from human TGFR1 around the phosphorylation site of Thr204.
-----------	---

Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

[Anti-Phospho-TGFBR1 \(T204\) Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.