

Anti-Phospho-DYRK1A/B (Y319/Y271) Antibody

Catalog Number: P00878

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Phospho-DYRK1A/B (Y319/Y271) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-DYRK1A/B (Y319/Y271) Antibody catalog # P00878. Tested in IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt: -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q13627/Q9Y463

Technical Details

Immunogen	A synthesized peptide derived from human DYRK1A/B around the phosphorylation site of Tyr319/271.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-40

Anti-Phospho-DYRK1A/B (Y319/Y271) Antibody

For Research Use Only. Not for use in diagnostic procedures.