

Anti-Phospho-HDAC4 (S632) Antibody

Catalog Number: P00971

About HDAC4

Histone deacetylase 4, also known as HDAC4, is a protein that in humans is encoded by the HDAC4 gene. It is mapped to 2q37.3. Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class II of the histone deacetylase/acuc/apha family. It possesses histone deacetylase activity and represses transcription when tethered to a promoter. This protein does not bind DNA directly, but through transcription factors MEF2C and MEF2D. It seems to interact in a multiprotein complex with RbAp48 and HDAC3.

Overview

Product Name	Anti-Phospho-HDAC4 (S632) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-HDAC4 (S632) Antibody catalog # P00971. Tested in WB, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P56524

Technical Details

Immunogen	A synthesized peptide derived from human HDAC4 around the phosphorylation site of Ser632.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-HDAC4 (S632) Antibody

For Research Use Only. Not for use in diagnostic procedures.