

Anti-Phospho-POLR2A (S5) Rabbit Monoclonal Antibody

Catalog Number: P01029-2


Overview

Product Name	Anti-Phospho-POLR2A (S5) Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-POLR2A (S5) Rabbit Monoclonal Antibody catalog # P01029-2. Tested in WB, IHC, ICC/IF, IP, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Monoclonal ABCE-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P24928


Technical Details

Immunogen	A synthesized peptide derived from human POLR2A
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 IP 1:20 FC 1:40


Anti-Phospho-POLR2A (S5) Rabbit Monoclonal Antibody (P01029-2) Images



Western blot analysis of Phospho-POLR2A (S5) expression in (1) HeLa cell lysate treated with Lambda Phosphatase; (2) HeLa cell lysate.



Immunohistochemical analysis of paraffin-embedded human spleen, using Phospho-POLR2A (S5) Antibody.



Immunofluorescent analysis of HeLa cells, using Phospho-POLR2A (S5) Antibody .

1 Publications Citing This Product

1. PubMed ID: 28848320, Heparanase mediates vascular endothelial growth factor gene transcription in high-glucose human retinal microvascular endothelial cells.

Visit bosterbio.com/anti-phospho-polr2a-s5-rabbit-monoclonal-antibody-p01029-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-POLR2A (S5) Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.