

Anti-Phospho-RPS6KA1/2/3/4 (S221/S227/S218/S232) Antibody

Catalog Number: P01058-4

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Phospho-RPS6KA1/2/3/4 (S221/S227/S218/S232) Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-Phospho-RPS6KA1/2/3/4 (S221/S227/S218/S232) Antibody catalog # P01058-4. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat, Monkey.
Application	IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q15418/P51812/Q15349/Q9UK32

Technical Details

Immunogen	A synthesized peptide derived from human RSK1/2/3/4 around the phosphorylation site of Ser221/227/s218/232.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

[Anti-Phospho-RPS6KA1/2/3/4 \(S221/S227/S218/S232\) Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.