

## Anti-Phospho-gp130/IL6ST (S782) Antibody

Catalog Number: P01216

### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Phospho-gp130/IL6ST (S782) Antibody                                                                                                                                                                                                                                                                                                                                                                                   |
| Reactive Species     | Human, Mouse                                                                                                                                                                                                                                                                                                                                                                                                               |
| Description          | Boster Bio Anti-Phospho-gp130/IL6ST (S782) Antibody catalog # P01216. Tested in WB, IHC applications. This antibody reacts with Human, Mouse.                                                                                                                                                                                                                                                                              |
| Application          | IHC, WB                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                       |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Uniprot ID           | P40189                                                                                                                                                                                                                                                                                                                                                                                                                     |

### Technical Details

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Immunogen           | A synthesized peptide derived from human CD130/gp130 around the phosphorylation site of Ser782. |
| Form                | Liquid                                                                                          |
| Concentration       | 500 ug/ml                                                                                       |
| Purification        | Protein A affinity purified.                                                                    |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-200                                      |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-gp130/IL6ST (S782) Antibody

For Research Use Only. Not for use in diagnostic procedures.