

Anti-Phospho-AMPK alpha (T172) PRKAA2 Rabbit Monoclonal Antibody

Catalog Number: P01420-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Phospho-AMPK alpha (T172) PRKAA2 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-AMPK alpha (T172) PRKAA2 Rabbit Monoclonal Antibody catalog # P01420-2. Tested in WB application. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal EHA-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13131/P54646


Technical Details

Immunogen	A synthesized peptide derived from human AMPK alpha
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography


Suggested Dilutions

WB 1:500-1:2000

Anti-Phospho-AMPK alpha (T172) PRKAA2 Rabbit Monoclonal Antibody (P01420-2) Images



All lanes use the Antibody at 1:4K dilution for 1 hour at room temperature.



Western blot analysis of Phospho-AMPK alpha (T172) expression in Mouse heart lysate.

[Anti-Phospho-AMPK alpha \(T172\) PRKAA2 Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.