

## Anti-Phospho-CDC25A (S124) Antibody

Catalog Number: P01433

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-Phospho-CDC25A (S124) Antibody  |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-Phospho-CDC25A (S124) Antibody catalog # P01433. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.   |
| Application          | IF, IHC, ICC, WB   |
| Clonality            | Polyclonal   |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.   |
| Host                 | Rabbit   |
| Uniprot ID           | P30304   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | A synthesized peptide derived from human CDC25A around the phosphorylation site of Ser124.                     |
| Form                | Liquid   |
| Concentration       | 500 ug/ml  |
| Purification        | Protein A affinity purified.   |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-200<br>Immunocytochemistry/Immunofluorescence, 1:50-200 |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



**Anti-Phospho-CDC25A (S124) Antibody**

For Research Use Only. Not for use in diagnostic procedures.