

## Anti-Phospho-RAD52 (Y104) Antibody

Catalog Number: P01580

### Overview

Product Name	Anti-Phospho-RAD52 (Y104) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho-RAD52 (Y104) Antibody catalog # P01580. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P43351

### Technical Details

Immunogen	A synthesized peptide derived from human RAD52 around the phosphorylation site of Tyr104.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-RAD52 (Y104) Antibody

For Research Use Only. Not for use in diagnostic procedures.