

## Anti-Phospho-PKC zeta (T560) PRKCZ Rabbit Monoclonal Antibody

Catalog Number: P01796

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Phospho-PKC zeta (T560) PRKCZ Rabbit Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                      |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Description          | Boster Bio Anti-Phospho-PKC zeta (T560) PRKCZ Rabbit Monoclonal Antibody catalog # P01796. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                            |
| Application          | IHC, WB                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Clonality            | Monoclonal GED-16                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | Q05513                                                                                                                                                                                                                                                                                                                                                                                                                                             |


### Technical Details

|               |                                                   |
|---------------|---------------------------------------------------|
| Immunogen     | A synthesized peptide derived from human PKC zeta |
| Isotype       | Rabbit IgG                                        |
| Form          | Liquid                                            |
| Concentration | 0.5mg/ml                                          |
| Purification  | Affinity-chromatography                           |


Suggested Dilutions

WB 1:500-2000  
IHC 1:50-200


## Anti-Phospho-PKC zeta (T560) PRKCZ Rabbit Monoclonal Antibody (P01796) Images



Western blot analysis of Phospho-PKC zeta (T560) expression in HeLa cell treated with Calyculin A lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

[Anti-Phospho-PKC zeta \(T560\) PRKCZ Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.