

Anti-Phospho Acetyl CoA Carboxylase (S79) (ACC1) ACACA Rabbit Monoclonal Antibody, Clone#RM270

Catalog Number: P01802


Overview

Product Name	Anti-Phospho Acetyl CoA Carboxylase (S79) (ACC1) ACACA Rabbit Monoclonal Antibody, Clone#RM270
Reactive Species	Human
Description	Boster Bio Anti-Phospho Acetyl CoA Carboxylase (S79) (ACC1) ACACA Rabbit Monoclonal Antibody, Clone#RM270 (Catalog # P01802). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM270
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13085


Technical Details

Immunogen	A phospho-peptide corresponding to human Phospho- Acetyl CoA Carboxylase (Ser79)
Cross Reactivity	This antibody reacts to human Acetyl CoA Carboxylase (ACC1) only when phosphorylated at Ser79. There is no cross-reactivity to Acetyl CoA Carboxylase without phosphorylation at Ser79.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:1000-1:2000 dilution WB: 1:1000-1:5000 dilution.

Anti-Phospho Acetyl CoA Carboxylase (S79) (ACC1) ACACA Rabbit Monoclonal Antibody, Clone#RM270 (P01802) Images



Western Blotting result Western Blot of A431 cell lysates, nontreated (-) or treated (+) with Lambda Protein Phosphatase (PP), using Anti-Phospho-Acetyl CoA Carboxylase (Ser79) Rabbit Monoclonal Antibody (Clone RM269) at a 1:1000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human breast cancer tissue sections using Anti-Phospho-Acetyl CoA Carboxylase (Ser79) Rabbit Monoclonal Antibody (clone RM270) at a 1:5000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho Acetyl CoA Carboxylase (S79) (ACC1) ACACA Rabbit Monoclonal Antibody, Clone#RM270

For Research Use Only. Not for use in diagnostic procedures.