

Anti-Phospho-PDHA1/2 (S293/S291) Antibody

Catalog Number: P01906-1

Overview

Product Name	Anti-Phospho-PDHA1/2 (S293/S291) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-PDHA1/2 (S293/S291) Antibody catalog # P01906-1. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P08559/P29803

Technical Details

Immunogen	A synthesized peptide derived from human PDHA1/2 around the phosphorylation site of Ser293/291.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200 Immunocytochemistry/Immunofluorescence, 1:50-200

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Phospho-PDHA1/2 (S293/S291) Antibody

For Research Use Only. Not for use in diagnostic procedures.