

Anti-Phospho-IKK Alpha/Beta (S176/S177) Antibody

Catalog Number: P01918-1

About CHUK/IKBKB

IKBKA(INHIBITOR OF KAPPA LIGHT POLYPEPTIDE GENE ENHANCER IN B CELLS, KINASE OF, ALPHA), also called CHUK, NFKBIKA, IKKA, IKK1, I-KAPPA-B, is a protein kinase that in humans is encoded by the CHUK gene. And it is a member of the protein kinase superfamily. And the protein contains a serine/threonine kinase domain, due to its ubiquitous expression in a broad array of tissues and high degree of conservation across species. By FISH, the IKBKA gene is mapped to chromosome 10q24. Phosphorylation of serine residues on the I-kappa-B proteins by kinases(IKBKA or IKBKB) marks them for destruction via the ubiquitination pathway, thereby allowing activation of the NF-kappa-B complex. IKK-alpha is part of the I kappa B kinase complex that plays an important role in regulating the NF-kappa B transcription factor.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Phospho-IKK Alpha/Beta (S176/S177) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-IKK Alpha/Beta (S176/S177) Antibody catalog # P01918-1. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O15111/O14920

Technical Details

--	--

Immunogen	A synthesized peptide derived from human IKK alpha/ beta around the phosphorylation site of Ser176.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200 Immunocytochemistry/Immunofluorescence, 1:50-200

[Anti-Phospho-IKK Alpha/Beta \(S176/S177\) Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.