

Anti-Phospho-Myosin Light Chain 2/MLC-2V/MYL2 (S19) Antibody

Catalog Number: P03215-1

Overview

Product Name	Anti-Phospho-Myosin Light Chain 2/MLC-2V/MYL2 (S19) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-Myosin Light Chain 2/MLC-2V/MYL2 (S19) Antibody catalog # P03215-1. Tested in IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P10916

Technical Details

Immunogen	A synthesized peptide derived from human MLC2 around the phosphorylation site of Ser19.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-Myosin Light Chain 2/MLC-2V/MYL2 (S19) Antibody

For Research Use Only. Not for use in diagnostic procedures.