

Anti-Phospho OSR1 (Thr185) Antibody

Catalog Number: P05141

Overview


Product Name	Anti-Phospho OSR1 (Thr185) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho OSR1 (Thr185) Antibody (Catalog# P05141). Tested in WB, IHC-P, IF (paraffin), ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC-P, WB, IF (paraffin)
Clonality	Polyclonal CA-20
Formulation	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95747

Technical Details


Immunogen	A peptide derived from human OSR1 around the phosphorylation site of Thr185. Immunogen sequence location: 151-200
Predicted Reactive Species	Human
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>Western Blot, 1:500 - 1:2000</p> <p>Immunohistochemistry, 1:100 - 1:300</p> <p>ELISA, 1:10000</p>

	Not yet tested in other applications.
--	---------------------------------------

Anti-Phospho OSR1 (Thr185) Antibody (P05141) Images



Western blot analysis of lysates from HepG2 cells treated with serum 20% 15', using OSR1 (Phospho-Thr185) Antibody. The lane on the right is blocked with the phospho peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using OSR1 (Phospho-Thr185) Antibody. The picture on the right is blocked with the phospho peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho OSR1 (Thr185) Antibody