

Anti-Phospho-RPS6KB2 (S423) Antibody

Catalog Number: P05845

Overview

Product Name	Anti-Phospho-RPS6KB2 (S423) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho-RPS6KB2 (S423) Antibody catalog # P05845. Tested in WB, ICC, IF applications. This antibody reacts with Human, Mouse.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UBS0

Technical Details

Immunogen	A synthesized peptide derived from human p70 S6 Kinase beta around the phosphorylation site of Ser423.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-RPS6KB2 (S423) Antibody

For Research Use Only. Not for use in diagnostic procedures.