

Anti-Phospho-Histone H3.3 (S31) Rabbit Monoclonal Antibody

Catalog Number: P06819-5

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Phospho-Histone H3.3 (S31) Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-Histone H3.3 (S31) Rabbit Monoclonal Antibody catalog # P06819-5. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal 28H13
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P84243


Technical Details

Immunogen	A synthesized peptide derived from human Phospho-Histone H3.3 (S31)
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions

WB 1:500-2000
IHC 1:50-200
ICC/IF 1:50-200
FC 1:50

Anti-Phospho-Histone H3.3 (S31) Rabbit Monoclonal Antibody (P06819-5) Images



Western blot analysis of Phospho-Histone H3.3 (S31) expression in HeLa cell treated with Calyculin A lysate.

Anti-Phospho-Histone H3.3 (S31) Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.