

## Anti-Phospho-IkB Epsilon/NFKBIE (S161) Antibody

Catalog Number: P07073

### About NFKBIE

Nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, epsilon, also known as NFKBIE, is a protein which in humans is encoded by the NFKBIE gene. Predicted to be involved in D-serine transport. Predicted to act upstream of or within negative regulation of transcription, DNA-templated. Predicted to be located in several cellular components, including Golgi apparatus; nuclear lumen; and perinuclear region of cytoplasm. Is expressed in several structures, including central nervous system; dorsal root ganglion; eye; genitourinary system; and thymus. Orthologous to human NFKBIE (NFKB inhibitor epsilon).

### Overview

Product Name	Anti-Phospho-IkB Epsilon/NFKBIE (S161) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-IkB Epsilon/NFKBIE (S161) Antibody catalog # P07073. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O00221

### Technical Details

Immunogen	A synthesized peptide derived from human I kappaB- epsilon around the phosphorylation site of Ser161.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

Immunocytochemistry/Immunofluorescence, 1:50-200

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-IkB Epsilon/NFKBIE (S161) Antibody

For Research Use Only. Not for use in diagnostic procedures.