

Anti-Phospho-FUNDC1 (S17) Antibody

Catalog Number: P08688

About FUNDC1

This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known.

Overview

Product Name	Anti-Phospho-FUNDC1 (S17) Antibody
Reactive Species	Human
Description	Boster Bio Anti-Phospho-FUNDC1 (S17) Antibody catalog # P08688. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8IVP5

Technical Details

Immunogen	A synthesized peptide derived from human FUNDC1 around the phosphorylation site of Serine 17.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-FUNDC1 (S17) Antibody

For Research Use Only. Not for use in diagnostic procedures.