

Anti-C-C chemokine receptor type 5 CCR5 Antibody Picoband®

Catalog Number: PA1016-2

About Ccr5

Human CC chemokine receptor 5 (CCR5), which is 17.5 kb from the CMKBR2 gene, encodes a 352-amino acid protein with a calculated molecular mass of 40,600 Da. It is organized into four exons and two introns. Human CCR5 is a co-receptor for macrophage-tropic strains of human immunodeficiency virus (HIV)-1 and is expressed by bone-marrow-derived cells. CCR5 is expressed on neurons, astrocytes, and microglia in the central nervous system, and on epithelium, endothelium, vascular smooth muscle and fibroblasts in other tissues. Functionally, CCR5 mediates the activation of cells by the chemokines macrophage inflammatory protein-1alpha, macrophage inflammatory protein-1beta, and RANTES, and serves as a fusion cofactor for macrophage-tropic strains of human immunodeficiency virus type 1.

Overview


Product Name	Anti-C-C chemokine receptor type 5 CCR5 Antibody Picoband®
Reactive Species	Mouse
Description	Boster Bio Anti-C-C chemokine receptor type 5 CCR5 Antibody catalog # PA1016-2. Tested in WB applications. This antibody reacts with Mouse. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Thimerosal, 0.05mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51682

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of mouse CCR5.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.

Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Mouse

Anti-C-C chemokine receptor type 5 CCR5 Antibody Picoband® (PA1016-2) Images



Anti-CCR5 antibody, PA1016-2, Western blotting
Lane 1: Mouse Lung Tissue Lysate
Lane 2: Mouse Intestine Tissue Lysate
Lane 3: Mouse Kidney Tissue Lysate

2 Publications Citing This Product

1. PubMed ID: 23063706, Li G, Xu H, Zhu S, Xu W, Qin S, Liu S, Tu G, Peng H, Qiu S, Yu S, Zhu Q, Fan B, Zheng C, Li G, Liang S. Brain Res Bull. 2013 Jan;90:79-87. Doi: 10.1016/J.Brainresbull.2012.10.002. Epub 2012 Oct 9. Effects Of Neferine On Ccl5 And Ccr5 Expression In...
2. PubMed ID: 22238429, Study on CCR5 analogs and affinity peptides

Visit bosterbio.com/anti-ccr5-antibody-pa1016-2-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C-C chemokine receptor type 5 CCR5 Antibody

For Research Use Only. Not for use in diagnostic procedures.