

Anti-Integrin alpha 1/ITGA1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA1045-1-Fluoro488

About ITGA1

Integrin alpha 1 (ITGA1) chain associates with the beta 1 (ITGB1) chain to form a heterodimer that functions as a dual laminin/collagen receptor in neural cells and hematopoietic cells. ITGA1 has a 206-amino acid I domain in its N-terminal half, followed by 3 divalent cation-binding sites and a C-terminal transmembrane domain with a short cytoplasmic tail. It also has 28 potential N-glycosylation sites. Human ITGA1 was expressed in a mouse fibroblast cell line as a 180-kD protein. ITGA1 is involved in the early remodeling of osteoarthritic cartilage and plays an essential role in the regulation of mesenchymal stem cell proliferation and cartilage production. It also plays an essential role in the regulation of MSC proliferation and cartilage production.

Overview

| | |
|----------------------|---|
| Product Name | Anti-Integrin alpha 1/ITGA1 Antibody Picoband® Fluoro488 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P56199 |

Technical Details

| | |
|------------------|--|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human Integrin alpha 1, different from the related mouse and rat sequences by one amino acid. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm |

Suggested Dilutions

Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 24643246, A Small Physiological Electric Field Mediated Responses of Extravillous Trophoblasts Derived from HTR8/SVneo Cells: Involvement of Activation of Focal Adhesion Kinase Signaling

Visit bosterbio.com/anti-integrin-alpha-1-antibody-pa1045-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin alpha 1/ITGA1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.