

Anti-Myelin Basic Protein/MBP Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PA1050-Fluoro488

About MBP

Myelin basic protein (MBP) is a major constituent of the myelin sheath of oligodendrocytes and Schwann cells in the central nervous system and the peripheral nervous system, respectively. It is most abundant in hemopoietic system and contains seven exons distributed over 32-34 kb. MBP isolated from MS brain may differ in charge microheterogeneity which would affect antigenic determinants. MBP is mapped to chromosome 18q22-23. Failure in this gene expression would be correlated in the central white matter with extrapyramidal system degeneration signs. Moreover, it is a candidate autoantigen in the disease multiple sclerosis.

Overview

Product Name	Anti-Myelin Basic Protein/MBP Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Pig, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P02686

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Myelin Basic Protein, identical to the related rat and mouse sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

17 Publications Citing This Product

1. PubMed ID: 10.1016/j.celrep.2015.01.037, Focal MMP-2 and MMP-9 Activity at the Blood-Brain Barrier Promotes Chemokine-Induced Leukocyte Migration

2. PubMed ID: 10.1016/j.brainresbull.2020.10.015, [Met5]-enkephalin preserves diffusion metrics in EAE mice

3. PubMed ID: 32515838, Meng FW, Jing XN, Song GH, Jie LL, Shen FF. Prox1 induces new lymphatic vessel formation and promotes nerve reconstruction in a mouse model of sciatic nerve crush injury. *J Anat.* 2020 Nov; 237(5):933-940. doi: 10.1111/joa.13247. Epub 2020 Jun 9. PMID: 32515838; PMCID:

Visit bosterbio.com/anti-myelin-basic-protein-antibody-pa1050-boster.html to see all 17 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Myelin Basic Protein/MBP Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.