

## Anti-NMDAR2A/GRIN2A Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PA1058-1-Fluoro594

### About GRIN2A

GRIN2A is also known as N-methyl-D-aspartate receptor channel, subunit epsilon-1 (NMDAR2A). This gene encodes a member of the glutamate-gated ion channel protein family. The encoded protein is an N-methyl-D-aspartate (NMDA) receptor subunit. NMDA receptors are both ligand-gated and voltage-dependent, and are involved in long-term potentiation, an activity-dependent increase in the efficiency of synaptic transmission thought to underlie certain kinds of memory and learning. These receptors are permeable to calcium ions, and activation results in a calcium influx into post-synaptic cells, which results in the activation of several signaling cascades. Disruption of this gene is associated with focal epilepsy and speech disorder with or without mental retardation. Alternative splicing results in multiple transcript variants.

### Overview

Product Name	Anti-NMDAR2A/GRIN2A Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q12879

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human NMDAR2A, different from the related mouse sequence by three amino acids, and from the related rat sequence by four amino acids.
Cross Reactivity	No cross reactivity with other proteins
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

1. PubMed ID: 10.1016/j.neuro.2016.01.005, Impaired hippocampal synaptic plasticity and NR2A/2B expression ratio in remifentanyl withdrawal rats

2. PubMed ID: 33099751, Liang Y, Ma Y, Wang J, Nie L, Hou X, Wu W, Zhang X, Tian Y. Leptin Contributes to Neuropathic Pain via Extrasynaptic NMDAR-nNOS Activation. Mol Neurobiol. 2020 Oct 25. doi:10.1007/s12035-020-02180-1. Epub ahead of print. PMID:33099751.

Visit [bosterbio.com/anti-nmdar2a-antibody-pa1058-1-boster.html](https://bosterbio.com/anti-nmdar2a-antibody-pa1058-1-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NMDAR2A/GRIN2A Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.