

## Anti-PP2A-alpha/PPP2CA Antibody Picoband® APC Conjugated

Catalog Number: PA1068-APC

### About PPP2CA

The catalytic subunit of human protein phosphatase 2A (PPP2CA) encodes a 309-amino acid polypeptide. It is localized to chromosome 5. The gene (approximately 30 kbp) is composed of seven exons and six introns. It is predicted to be important for phosphatase enzymatic activity. Methylation of the C-terminal leucine residue (Leu309) of protein serine/threonine phosphatase 2A catalytic subunit (PP2AC) is known to regulate catalytic activity in vitro. Furthermore, PP2A has a fundamental role in cardiac function, and suggests that disturbances in protein phosphatase expression and activity may cause or exacerbate the course of cardiac diseases.

### Overview

Product Name	Anti-PP2A-alpha/PPP2CA Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P67775

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PP2A-alpha, identical to the related rat and mouse sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 19891536, Tian Hp, Qiu Tz, Zhao J, Li Lx, Guo J. Brain Inj. 2009 Dec;23(13-14):1073-80. Doi: 10.3109/02699050903379388. Sphingomyelinase-Induced Ceramide Production Stimulate Calcium-Independent Jnk And Pp2A Activation Following Cerebral Ischemia.

Visit [bosterbio.com/anti-pp2a-alpha-antibody-pa1068-boster.html](http://bosterbio.com/anti-pp2a-alpha-antibody-pa1068-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PP2A-alpha/PPP2CA Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.