

## Anti-IL8/CXCL8 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PA1353-Fluoro550

### About CXCL8

Interleukin-8 (IL-8) is a chemokine produced by macrophages and other cell types such as epithelial cells. It is also synthesized by endothelial cells, which store IL-8 in their storage vesicles, the Weibel-Palade bodies. In humans, the interleukin-8 protein is encoded by the IL8 gene. Interleukin-8 (IL8) is a member of the CXC chemokine family (Hull et al., 2001). IL-8 is believed to play a role in the pathogenesis of bronchiolitis, a common respiratory tract disease caused by viral infection. This gene and other ten members of the CXC chemokine gene family form a chemokine gene cluster in a region mapped to chromosome 4q. The genes for IL8 have been co-localized on a 335-kb genomic fragment.

### Overview

Product Name	Anti-IL8/CXCL8 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P10145

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human IL8.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 46 Publications Citing This Product

---

1. PubMed ID: 26949603, Inhibition of zymosan-induced cytokine and chemokine expression in human corneal fibroblasts by triptolide
2. PubMed ID: 28726780, FOXQ1 regulates senescence-associated inflammation via activation of SIRT1 expression
3. PubMed ID: 25317156, Neuroprotective effect of pretreatment with ganoderma lucidum in cerebral ischemia/reperfusion injury in rat hippocampus

Visit [bosterbio.com/anti-il-8-antibody-pa1353-boster.html](http://bosterbio.com/anti-il-8-antibody-pa1353-boster.html) to see all 46 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL8/CXCL8 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.