

## Anti-Interleukin-10 IL10 Antibody Picoband®

Catalog Number: PA1354

### About IL10

Interleukin-10 (IL-10 or IL10), also known as human cytokine synthesis inhibitory factor (CSIF), is an anti-inflammatory cytokine. In humans IL-10 is encoded by the IL10 gene. It is capable of inhibiting synthesis of pro-inflammatory cytokines like IFN-gamma, IL-2, IL-3, TNFalpha and GM-CSF made by cells such as macrophages and regulatory T-cells. IL-10 also displays potent abilities to suppress the antigen presentation capacity of antigen presenting cells. Kim et al. (1992) showed that the mouse Il 10 gene contains 5 exons and spans about 5.2 kb of genomic DNA. Eskdale et al. (1997) mapped the IL10 gene to the junction between 1q31 and 1q32.

### Overview


Product Name	Anti-Interleukin-10 IL10 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-Interleukin-10 IL10 Antibody catalog # PA1354. Tested in WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg Thimerosal, 0.05mg NaN3. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22301

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human IL10.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human

## Anti-Interleukin-10 IL10 Antibody Picoband® (PA1354) Images



Anti-IL-10 antibody, PA1354, Western blotting  
Lane 1: Recombinant Human IL-10 Protein 10ng  
Lane 2: Recombinant Human IL-10 Protein 5ng  
Lane 3: Recombinant Human IL-10 Protein 2.5ng

## 17 Publications Citing This Product

1. PubMed ID: 25245263, Zhan X, Jia L, Niu Y, Qi H, Chen X, Zhang Q, Zhang J, Wang Y, Dong L, Wang C. Biomaterials. 2014 Dec;35(38):10046-57. Doi: 10.1016/J.Biomaterials.2014.09.007. Epub 2014 Sep 22. Targeted Depletion Of Tumour-Associated Macrophages By An Alendronate-...
2. PubMed ID: 26083117, Protective Effect of the Total Saponins from Rosa laevigata Michx Fruit against Carbon Tetrachloride-Induced Liver Fibrosis in Rats
3. PubMed ID: 22594905, Zhong B, Shen H, Sun X, Wang H, Zhang Y, Sun Z. Cancer Biother Radiopharm. 2012 May;27(4):252-8. Doi: 10.1089/Cbr.2011.1105. Additive Effects Of Ulinastatin And Docetaxel On Growth Of Breast Cancer Xenograft In Nude Mice And Expression Of Pge2, Il...

Visit [bosterbio.com/anti-il-10-antibody-pa1354-boster.html](http://bosterbio.com/anti-il-10-antibody-pa1354-boster.html) to see all 17 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-10 IL10 Antibody

For Research Use Only. Not for use in diagnostic procedures.