

## Anti-Presenilin 2/PSEN2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PA1358-Fluoro594

### About PSEN2

Presenilin-2 is a protein that in humans is encoded by the PSEN2 gene. Kovacs et al. (1996) demonstrated that the expression patterns of PS1 and PS2 in the brain are extremely similar to each other and that messages for both are primarily detectable in neuronal populations. Immunochemical analyses indicated that PS1 and PS2 are similar in size and localize to similar intracellular compartments (endoplasmic reticulum and Golgi complex). Li et al. (1997) demonstrated that wildtype PS1 and PS2 are localized to the nuclear membrane, its associated interphase kinetochores, and the centrosomes. In melanocytic cells PSEN2 gene expression may be regulated by MITF.

### Overview

Product Name	Anti-Presenilin 2/PSEN2 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P49810

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Presenilin 2, different from the related mouse and rat sequences by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Presenilin 2/PSEN2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.