

## Anti-Insulin-like growth factor I IGF1 Antibody Picoband® APC Conjugated

Catalog Number: PA1374-APC

### About IGF1

Insulin-like growth factor 1 (IGF-1) also known as somatomedin C or mechano growth factor is a protein that in humans is encoded by the IGF1 gene. IGF-1 is a hormone similar in molecular structure to insulin. It plays an important role in childhood growth and continues to have anabolic effects in adults. A synthetic analog of IGF-1, mecasermin is used for the treatment of growth failure. IGF-1 consists of 70 amino acids in a single chain with three intramolecular disulfide bridges. IGF-1 has a molecular weight of 7649 daltons. Justice et al. (1990) placed the mouse IGF1 gene on chromosome 10.

### Overview

Product Name	Anti-Insulin-like growth factor I IGF1 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P05019

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human IGF-1, different from the mouse and rat sequences by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 10 Publications Citing This Product

1. PubMed ID: 33617969, Peng MF, Tian S, Song YG, Li CX, Miao MS, Zhen-Ren, Li M. Effects of total flavonoids from *Eucommia ulmoides* Oliv. leaves on polycystic ovary syndrome with insulin resistance model rats induced by letrozole combined with a high-fat diet. *J Ethnopharmacol.* 2021 Feb
2. PubMed ID: 27779666, Effects of the IGF-1/PTEN/Akt/FoxO signaling pathway on male reproduction in rats subjected to water immersion and restraint stress
3. PubMed ID: 27833113, Elucidating a molecular mechanism that the deterioration of porcine meat quality responds to increased cortisol based on transcriptome sequencing

Visit [bosterbio.com/anti-igf-1-antibody-pa1374-boster.html](https://bosterbio.com/anti-igf-1-antibody-pa1374-boster.html) to see all 10 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Insulin-like growth factor I IGF1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.