

## Anti-IL23 P19/IL23A Antibody Picoband®

Catalog Number: PA1470

### About IL23a

IL-23, Interleukin-23, also known as Interleukin-23 subunit alpha1n (IL23A), is a heterodimeric cytokine consisting of two subunits, one called p40, which is shared with another cytokine, IL-12, and another called p19 (the IL-23 alpha subunit). The International Radiation Hybrid Mapping Consortium mapped the IL-23 gene to chromosome 12. IL-23 is an important part of the inflammatory response against infection. It promotes upregulation of the matrix metalloproteinase MMP9, increases angiogenesis and reduces CD8+ T-cell infiltration.

### Overview


Product Name	Anti-IL23 P19/IL23A Antibody Picoband®
Reactive Species	Mouse
Description	Boster Bio Anti-IL23 P19/IL23A Antibody catalog # PA1470. Tested in WB applications. This antibody reacts with Mouse. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg Thimerosal, 0.05mg NaN <sub>3</sub> . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9EQ14

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of mouse IL23 P19, different from the related rat sequence by two amino acids.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Mouse

## Anti-IL23 P19/IL23A Antibody Picoband® (PA1470) Images



All lanes: Anti-IL23 P19 antibody, PA1470  
Lane 1: Recombinant Mouse IL-23 Protein 10ng  
Lane 2: Recombinant Mouse IL-23 Protein 5ng  
Lane 3: Recombinant Mouse IL-23 Protein 2.5ng

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL23 P19/IL23A Antibody

For Research Use Only. Not for use in diagnostic procedures.