

## Anti-STAT4 Antibody Picoband® HRP Conjugated

Catalog Number: PA1692-HRP

### About STAT4

STAT4 is a human gene. The protein encoded by this gene is a member of the STAT protein family of transcription factors. STAT4 is mapped to chromosome 2q32.2-q32.3 by fluorescence in situ hybridization. STAT4 is phosphorylated in response to interleukin-12 and is essential for IL12 signal transduction. Transactivation, electrophoretic mobility shift, and RNA interference analyses showed that Ikaros bound to the STAT4 promoter and was involved in regulation of STAT4 in human T cells.

### Overview

Product Name	Anti-STAT4 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q14765

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human STAT4, identical to the related rat and mouse sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-STAT4 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.