

## Anti-CYP11A1 Antibody Picoband® APC Conjugated

Catalog Number: PA1698-1-APC

### About Cyp11a1

CYP11A1 (Cytochrome p450, family 11, subfamily A, polypeptide 1), also called CYTOCHROME P450SCC, CYTOCHROME P450C11A1 or CYP11A, is a mitochondrial enzyme associated with the conversion of cholesterol to pregnenolone. CYP11A1 is a member of the cytochrome P450 superfamily of enzymes. The CYP11A1 gene is mapped on 15q24.1. Expression of the CYP11A1 gene may play a role in skin physiology and pathology and that cutaneous proopiomelanocortin activity may be autoregulated by a feedback mechanism involving glucocorticoids synthesized locally by Slominski et al. Using in vitro studies, CYP11A1 catalyzed the side-chain cleavage of 7-dehydrocholesterol to form 7-dehydropregnenolone. In addition, CYP11A1 catalyzed the metabolism of biologically inert vitamin D<sub>3</sub>, which is formed from 7-dehydrocholesterol, to form 2 hydroxylated products, 20-hydroxyvitamin D<sub>3</sub> and 20, 22-dihydroxyvitamin D<sub>3</sub>. Mutations in the CYP11A1 gene cause congenital adrenal insufficiency with partial or complete 46, XY sex reversal.

### Overview

Product Name	Anti-CYP11A1 Antibody Picoband® APC Conjugated
Reactive Species	Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9QZ82

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of mouse CYP11A1, different from the related rat sequence by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm

	Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 23555723, Abnormal Apoptosis of Trophoblastic Cells Is Related to the Up-Regulation of CYP11A Gene in Placenta of Preeclampsia Patients

Visit [bosterbio.com/anti-cyp11a1-antibody-pa1698-1-boster.html](http://bosterbio.com/anti-cyp11a1-antibody-pa1698-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CYP11A1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.